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## medic T02

AET-R161

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**BY Xblitz.**  
**neno**  
www.neno.pl  
inspired by Children  
designed by Parents.

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**Autoryzowany przedstawiciel:**  
Shanghai International Holding Corp. GmbH (Europe).  
Address: Eiffelstrasse 80, 20537 Hamburg, Germany

**Wytwórcą:**  
Alicn Medical Shenzhen, Inc  
4/F, B Building,  
Shenfabao Modern Optical Factory,  
Kengzi Street Pingshan District, 518122 Shenzhen City,  
People's Republic of China

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## medic T02 AET-R161

USER'S MANUAL

**Dear Customer,**  
thank you for purchasing NENO MEDIC T02 AET-R161. The device you received is a non-contact infrared, reusable thermometer. The device can be used both in home environment and by professionals. The thermometer measures body temperature of adults and children that are over 3 months old. **Before using the device, please refer to the instructions below. Keep this manual so you can refer to it later.**

### 1. INTRODUCTION

1. Try the device on yourself to understand how it works.
2. To avoid the temperature of the surroundings influencing the measurement, the thermometer should be kept at room temperature (16-35°C) for at least half an hour before the measurement is taken.
3. The person whose temperature is measured should stay inside a room with moderate temperature (16-35°C) in it for at least 20 minutes before taking the measurement.
4. After doing physical exercise, wait for 30 minutes before measuring the temperature.
5. Wipe the skin on the forehead that has its temperature measured by the thermometer and make sure that the hair does not obscure the path between the forehead and the sensor.
6. If you take multiple measurements in a row, wait 5 seconds between each measurement and take them from different distances to achieve the most accurate results.
7. Make sure that the sensor is clean before each measurement.
8. The standard body temperature, both for adults and for children, is about 36.6 degrees Celsius. Consult your doctor if the measured temperature is much higher or lower than that.

### 2. HOW TO USE

- The thermometer measures the temperature of a person's forehead and operates at a distance of 1.5 cm to 5 cm.
1. Turn on the thermometer by pressing the **ON/OFF** button. When you hear a beep and **---**°C is displayed on the screen, the unit is ready to use.
  2. Place the thermometer at an appropriate distance from the forehead of the person that is being measured, with the sensor pointing towards the forehead. Press **▶** and wait a second for a short beep. The temperature measurement will be shown on the display screen. If you want to perform multiple tests in a row, wait 5 seconds between measurements. This procedure can be repeated as many times as you want.
  3. To turn off the device, hold **▶** for about 5 seconds. The thermometer will turn off automatically after a minute of inactivity.
  4. If the measurement is below 32°C, the "Lo" message will appear on the display and the thermometer will sound two beeps.
  5. If the measurement is above 37.8°C, the thermometer will sound six beeps.
  6. If the measurement is above 42.2°C, the "Hi" message appears on the display and the thermometer sound two beeps.


### 3. CHANGE THE UNIT OF MEASUREMENT

1. To change the unit of measurement from Celsius to Fahrenheit, hold down the **▶** button for 8 seconds while the thermometer is off until the screen displays **---**°C. Press **▶** to change the measurement unit. The thermometer will turn off automatically after 4 seconds.

### 4. MEASUREMENT MEMORY


1. Make sure the thermometer is turned off.
2. Hold **▶** for 4 seconds to turn on memory mode. The last measurement made by the thermometer will be displayed.
3. Press the **▶** button to switch between measurements saved in device memory. The higher the number displayed before the temperature, the older the measurement displayed.
4. Thermometer can show 32 last measurements.
5. In memory mode, the thermometer turns off automatically after 12 seconds of inactivity.

### 5. DESCRIPTION OF THE DEVICE

### 6. DISPLAY SCREEN

1. **MEASURING POSITION** - indicates the location (patient's forehead) that should be aimed at with the thermometer.
2. **TEMPERATURE DATA** - shows temperature measurement taken.
3. **TEMPERATURE UNIT** - shows whether the temperature is measured in Celsius or Fahrenheit degrees.
4. **BATTERY SYMBOL** - Appears when batteries need to be replaced with new ones.
5. **MEMORY SYMBOL** - shows that the measurement currently being viewed is an archival measurement stored in the device's memory.

### 7. BATTERY REPLACEMENT

### 8. SAFETY MEASURES

1. Do not use the thermometer in rooms where the temperature is below or above 35°C or where the humidity is below 15% RH or above 80%RH. These conditions may cause measurement errors.
2. Do not expose the thermometer to temperatures below -20°C or above 65°C or humidity that is below 15%RH or above 95%RH. These conditions may damage the device.
3. If the device sensor is damaged or loose, contact service immediately.
4. Do not touch or blow on the sensor. This may lead to incorrect measurements.
5. Do not make a self-diagnosis with the resulting measurements. If the results are of concern to you, consult a doctor.
6. If the sensor is dirty, gently wipe it with soft material. Cleaning the sensor with paper towels or toilet paper can scratch and damage it.
7. Incorrect installation of the battery in the device will lead to excessive battery heating and damage to the device.
8. Do not expose the device to liquids. The device is not waterproof.
9. The sensor is made of a special type of glass, which is difficult to damage. If you notice any damage, stop using the device immediately and contact the retailer.

10. Do not use your mobile phone near the thermometer that is turned on. Portable communication equipment can affect the measurement of medical electronic equipment.
11. Do not modify the device yourself and do not disassemble it. This will result in the loss of the warranty.
12. Keep the thermometer away from children. Some parts of the device are so small that the child can swallow them. If the child accidentally swallows the battery or protective film of the device, consult a doctor immediately.
13. If the device has been modified, appropriate tests must be carried out to ensure that the device is functioning properly. If you experience problems, please consult the service immediately.
14. If there is a temperature difference between the room where the thermometer is stored and the room where the measurement is made, make sure that the thermometer is kept for at least 30 minutes in the measuring room before the first measurement is carried out. This will help avoid irregularities.
15. The device does not require periodic calibration. In case of measurement errors, please contact the seller.
16. The person using the thermometer should not touch the thermometer and the patient at the same time. This can lead to measurement errors.
17. Do not maintain or clean the device while it is in use.
18. The patient can measure and replace the batteries by themselves by following the instructions in this manual.
19. The patient can maintain the device on their own.

### 9. ERROR MESSAGES

message	problem	solution
Hi	Temperature above 42.2°C	Use the thermometer only under the conditions specified in the manual. If the error persists, consult the service point.
Lo	Temperature below 32°C	Use the thermometer only under the conditions specified in the manual. If the error persists, consult the service point.
Err	Sensor not attached	Contact service
ErE	EEPROM memory isn't functioning properly	Contact service
	Discharged batteries	Replace the batteries with a new set
ErH	Ambient temperature too high	Go to a room with a temperature of 16-35°C
ErL	Ambient temperature too low	Go to a room with a temperature of 16-35°C

### 10. MAINTENANCE, STORAGE AND CALIBRATION OF THE DEVICE

- Maintenance:**
1. Remove stains from the device using soft, dry cloth.
  2. Gently wipe off the stains on the device with a delicate material, use a cotton swab or cloth with alcohol to gently wipe the sensor and display. After use, the product can be sterilized with a piece of cotton or a delicate cloth soaked in medical alcohol (75% alcohol content).
  3. Do not wash the appliance with water or cleaning agents containing corrosive agents or benzene.
  4. Do not submerge in any liquids.
- Storage:**
1. The device must not be stored in places where: the temperature is too high or low, the moisture is too high or too low, the device is exposed directly to sunlight, near the flow of electric current or in dusty rooms. These conditions can damage the device or cause inaccurate measurements.
  2. Remove the batteries from the device if it will not be used for a long time.
- Calibration:**
1. The thermometer is calibrated during the manufacturing process. If the thermometer is used in accordance with the instructions, it does not require periodic calibration. If you doubt the accuracy of the measurement, please contact the retailer that sold you the device or service center.
  2. Do not attempt to modify or disassemble the thermometer yourself.

### 11. SPECIFICATION

Measurement mode: Forehead  
 Measurement units: Celsius (°C) or Fahrenheit (°F)  
 Operating temperature: 16 to 35°C  
 Operating humidity: 15 to 80%  
 Storage temperature: -20°C to 55°C  
 Storage humidity: 15 to 93%  
 Measurement distance: 15 to 50mm from the measuring point  
 Measurement range : 32 to 42.2°C  
 Measurement accuracy : +/- 0.2° (in 35-42°C);  
 +/- 0.3°C in other ranges  
 Memory: 32 last measurements  
 Dimensions: 150x37x16,3mm  
 Weight: 51g (without batteries)  
 Batteries: 2xAAA (DC 3V)  
 Auto power off : After 1 minute of inactivity  
 Product Life: 5 years  
 Measure site: Forehead  
 Reference site: Armpit

### 12. EXPLANATION OF SYMBOLS

The following symbols may appear in the manual, on the thermometer, or on the rest of the set.

	Marks information that is mandatory to learn before you start using your device.
	Follow the instructions.
	BF-type part application
	<b>Warning:</b> Read the included documentation.
	Do not dispose of the product in a mixed municipal waste container. Dispose of the product in accordance with the disposal guidelines for electronic devices of this type.
	Transport packaging should be kept away from rain.
	Transport packaging should not be exposed to sunlight.
	Indicates the side of the transport packaging that should be facing upwards.
	The contents of the transport packaging are delicate, the packaging should be treated with extreme care.
	Indicates the temperature range at which the transport packaging should be stored.
	Non-ionizing electromagnetic radiation
	Don't roll.
	LOT Number
	Production date
	The product cannot be used after the date listed next to this icon.
	The product is resistant to solid objects with a diameter of 12.5mm or more and against vertically falling water droplets when the device is raised at an angle of up to 15 degrees.
	Manufacturer
	CE mark: The product complies with the Medical Device Directive 93/42/EEC
	Authorized representative in the European Community.

### 13. ELECTROMAGNETIC COMPATIBILITY

The ME DEVICE is adapted for use in home and similar environments for healthcare purposes.

**Warning:** Do not start the device near high-frequency surgical devices or in magnetic resonance rooms with RF shielding, as there are high amounts of electromagnetic disturbances in these locations.

**Note:** Using the device near other equipment may lead to measurement irregularities and should be avoided. If the use of two devices in each other's vicinity is unavoidable, both devices should be carefully monitored to ensure that they are working properly.

If any: a list of all cables and maximum cable length (if applicable), transducers and other ACCESSORIES which are replaceable by THE RESPONSIBLE ORGANIZATION and which may affect the compliance with the standards of ME EQUIPMENT or ME SYSTEM with the requirements set out in Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES can be specified in general (e.g. shielded cable) or specifically (e.g. specifying the MANUFACTURER and TYPE OF EQUIPMENT).

If any: The functioning of medical equipment or medical system has been defined as ESSENTIAL PERFORMANCE and a description the changes that the OPERATOR may experience if ESSENTIAL PERFORMANCE is disturbed by EM DISTURBANCES (the term ESSENTIAL PERFORMANCE does not need to be used).

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard for EM disturbances for the period of the device's expected service life.
2. Guidance and manufacturer's declaration - electromagnetic emissions and immunity.

Table 1

Guidance and manufacturer's declaration - electromagnetic emission	
Emission test	Standard met
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	N/A
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A

Table 2

Guidance and manufacturer's declaration - electromagnetic immunity		
Immunity test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/-2 kV +/-4 kV, +/-8kV, +/-15kV air	+/- 8 kV contact +/-2 kV +/-4 kV, +/-8kV, +/-15kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply cables: +/- 2kV Input/output cables: +/- 1kV	N/A
Surges IEC 61000-4-5	Line(s) to line(s): +/- 1kV Line(s) to earth: +/- 2kV Repetition frequency 100 kHz	N/A
Voltage dips, short interruptions and voltage variations in power supply cables IEC 61000-4-11	0% 0.5cycle At 0°, 45°, 90°, 135°, 180°, 225°, 275° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycles	N/A
Magnetic field frequency IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted radio frequencies IEC 61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	N/A
Radiated radio waves IEC 61000-4-3	10 V/m80MHz - 2,7GHz80% AM at 1kHz	10 V/m 80MHz - 2,7GHz 80% AM at 1kHz

NOTE: Ut is the AC mains voltage before applying the test level.

Table 3

Guidance and manufacturer's declaration - electromagnetic immunity						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Module (dB)	Distance (m)	Measured immunity level (dB)
385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	380-390	GMSR 460, FRS 460	FM +/- 5 kHz deviation 1 kHz sine	2	0.3	28
710 745 780	704-787	LTE band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
810 870 930	800-960	GSM 800/900, TETRA 800, iDEN 800, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
1720 1845 1970	1700-1990	GSM 1800, CDMA 1900, DECT, LTE Band 1, 3, 4, 25, UMTS	Pulse modulation 217Hz	2	0.3	28
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Pismo 7	Pulse modulation 217Hz	2	0.3	28
5240 5240 5785	5100-5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9

Radiated RF IEC 61000-4-3 (test specifications for ENCLOSED PORT (IMMUNITY) to RF wireless communication equipment).

### 14. WARRANTY CARD

The product is covered by a 24-month warranty. The terms of the guarantee can be found at: <https://neno.pl/gwarancja>

Complaints should be reported using the complaint form located at: <https://neno.pl/kontakt>

Specifications and contents of the kit are subject to change without notice. We apologise for any inconvenience.

KGK Trade declares that this device NENO Medic T02 is compliant with significant requirements of directive 2014/53/EU. The text of this declaration may be found on the website: <https://neno.pl/download/DOC/deklaracja-CE-Neno-Medic-T02.pdf>

The crossed out trash can symbol indicates that unusable electrical or electronic devices, its accessories (such as power supplies, cords) or components (for example batteries, if included) cannot be disposed of alongside with household waste. In order to dispose of the devices' its components (for example, batteries) deliver the device to the collection point, where it will be accepted free of charge. Disposal is subject to the recast version of the WEEE Directive (2012/19 / EU) and the Directive on batteries and accumulators (2006/66 / EC). Proper disposal of the device prevents degradation of the natural environment. Information about the collection points of the facilities is issued by the competent local authorities. Incorrect disposal of waste is subject to penalties provided for by the law in force in the given area.

**Importer:**  
KGK Trade sp. z o.o. sp. k.  
Os. Urocz 12  
31-953 Cracow  
Poland

**BY Xblitz.**  
**neno**  
www.neno.pl  
Inspired by Children,  
Designed by Parents

Made in P.R.C.

**EU Authorized Representative:**  
Shanghai International Holding Corp. GmbH (Europe).  
Address: Eiffelstrasse 80, 20537 Hamburg, Germany

**Manufacturer:**  
Alican Medical Shenzhen, Inc.  
4/F, B Building,  
ShenFubao Modern Optical Factory,  
Kengzi Street Pingshan District, 518122 Shenzhen City,  
People's Republic of China

The quality of the product is guaranteed by the TÜV SÜD Product Service GmbH

Last update of the user's manual: 30.07.2020